

How to submit samples or libraries?

→ complete our [submission form](https://ige3.lims.unige.ch/app/) (3 web pages) - <https://ige3.lims.unige.ch/app/>

All samples / libraries inside a submission must follow the same protocol.

If you have different types of submissions (e.g., RNA for RNA-seq and DNA for ChIP-seq), submit separate requests for each.

1 – « Contact informations » page

→ Enter the PI (Principal Investigator) and user information.

Note:

- Do not use special characters - including accented letters (e.g., é, ï).
- The .fastq files will be sent to the user.

2 – « Service » page

→ Specify the sequencing request (read length / depth) w/wo the analysis.

With:

Run type → select the desired read length for your submission.

Total reads needed (multiple of 200 mio reads +/- 20%)

→ this number applies to the entire submission

→ to calculate the correct value:

- multiply the number of samples by the reads required per sample.
- choose the closest multiple of 200 million reads.

Bioinformatics analysis → check the box if you want a bioinformatics analysis.

Type of request → indicate whether you are submitting **samples** (RNA, DNA, cells) or **pre-prepared libraries**.

3 – « Samples » page

→ Provide the characteristics common to all samples.

With:

Library preparation type → select the appropriate protocol to apply to your samples.

or « Libraries » page

→ Provide the characteristics common to all libraries.

With:

Library preparation kit chemistry → select the library prep kit chemistry used to construct the libraries.

→ Fill out the sample / library inventory table (direct entry or .csv file import – download the template).

With:

Tube label → ensure that the tube label matches what is physically written on the tube (numbers are preferred).

Sample name → the sample names you provide will be used to create the .fastq filename and can differ from tube label if needed. Do not use special characters and spaces (underscores are acceptable substitutes).

Note regarding the **index information**:

- Download the list of indexes included in our database (button "Indexes in the database")
- If your index is listed, only fill the *1st index name (i7)* column (and the *2nd index name (i5)* column if applicable).
- Else, provide the sequence in the *1st index sequence* column (and the *2nd index sequence* column if applicable). Be careful on the sequence orientation!

Then confirm your request:

→ Once the submission form is completed (the "SEND" button), you will receive an email which you must reply to in order to confirm your submission. If you do not receive an email, check your spam or contact us.

Available options for the Single Cell Sequencing:

	Submission form						
	« Service » page			« Samples » page		« Libraries » page	
You want to perform:	Run type	Total reads needed (multiple of 200 mio reads +/- 20%)		Sample type	Library preparation type	Library preparation kit chemistry	1st / 2nd index name
Chromium / Universal 3' Gene Expression							
a single-cell <u>whole transcriptome</u> profiling on <u>fresh or frozen</u> samples, w/wo on-chip <u>multiplexing</u> .	Paired End Read - 28 + 90 cycles	for GEX library for CSP library	20'000 read pairs/cell 5'000 read pairs/cell	specify if: single cells single nuclei	single-cell <u>RNA-seq</u> (10x genomics)	Chromium GEM-X 3' - v4 with multiplexing: Chromium GEM-X 3' 4-plex (OCM) - v4 previous kits: Chromium Next GEM 3' - v3.1	for GEX library SI-TT-** for CSP library SI-NT-**
Chromium / Flex Gene Expression							
a single-cell sequencing on <u>fixed</u> - even FFPE - samples (human, mouse).	Paired End Read - 28 + 90 cycles		10'000 read pairs/cell	<u>fixed</u> single cells	single-cell <u>RNA-seq Flex</u> (10x genomics)	Chromium GEM-X Flex Chromium GEM-X Flex – v2 (Apex) previous kits: Chromium Next GEM Flex	SI-TS-**
Chromium / Universal 5' Gene Expression							
a single-cell <u>gene expression</u> and <u>immunoprofiling</u> on <u>fresh or frozen</u> samples (human, mouse), w/wo on-chip <u>multiplexing</u> .	Paired End Read - 28 + 90 cycles	for GEX library for V(D)J library for CSP library for CRISPR library	20'000 read pairs/cell 5'000 read pairs/cell 5'000 read pairs/cell 5'000 read pairs/cell	specify if: single cells single nuclei	single-cell <u>RNA-seq & immunoprofiling</u> (10x genomics)	Chromium GEM-X 5' - v3 with multiplexing: Chromium GEM-X 5' 4-plex (OCM) - v3 previous kits: Chromium Next GEM 5' - v2	for GEX library SI-TT-** for V(D)J library SI-TT-** for CSP library SI-TN-** for CRISPR library SI-TT-**
Chromium / Epi Multiome ATAC + Gene Expression							
a single-cell <u>gene expression</u> and <u>epigenomic</u> profile from the same nucleus.	=> Proceed of two submissions! Add "_GEX" and "_ATAC" suffixes to the GEX and ATAC sample names respectively.						
submission 1:	Paired End Read - 28 + 90 cycles	for GEX library	20'000 read pairs/nucleus	single nuclei	single-cell <u>Multiome GEX</u> (10x genomics)	Chromium GEM-X Epi Multiome - GEX previous kits: Chromium Next GEM Multiome - GEX	for GEX library SI-TT-**
submission 2:	Paired End Read - 50 + 50 cycles	for ATAC library	25'000 read pairs/nucleus	single nuclei	single-cell <u>Multiome ATAC</u> (10x genomics)	Chromium GEM-X Epi Multiome - ATAC previous kits: Chromium Next GEM Multiome - ATAC	for ATAC library SI-NA-**
Chromium / Epi ATAC							
a single-cell <u>epigenomic</u> profile.	Paired End Read - 50 + 50 cycles		25'000 read pairs/nucleus	single nuclei	single-cell <u>ATAC</u> (10x genomics)	Chromium Next GEM ATAC - v2	SI-NA-**

GEX = Gene Expression
CSP = Cell Surface Protein
OCM = On-Chip Multiplexing

Available options for the Spatial Gene Expression:

Submission form						
« Service » page		« Samples » page		« Libraries » page		
You want to perform:	Run type	Total reads needed (multiple of 200 mio reads +/- 20%)	Sample type	Library preparation type	Library preparation kit chemistry	1st / 2nd index name
Visium / HD 3' Spatial Gene Expression a single-cell <u>whole transcriptome spatial</u> profiling, on <u>fresh or frozen</u> samples.	Paired End Read - 43 + 75 cycles (Visium 3')	550'000'000 reads pairs	single cells	Spatial Transcriptomics - <u>Visium HD 3'</u>	Spatial Transcriptomics - Visium HD 3'	SI-TT-**
Visium / HD Spatial Gene Expression a single-cell <u>whole transcriptome spatial</u> profiling, on <u>fresh, frozen or fixed</u> - even FFPE - samples (human, mouse).	Paired End Read - 43 + 50 cycles (Visium)	for fully-covered <u>6.5 mm</u> Capture Area: for FFPE sample 275'000'000 reads pairs for fixed sample 500'000'000 reads pairs for fresh/frozen sample 700'000'000 reads pairs for fully-covered <u>11 mm</u> Capture Area: for FFPE sample 825'000'000 reads pairs for fixed sample 1'500'000'000 reads pairs for fresh/frozen sample 2'100'000'000 reads pairs	specify if: single cells <u>fixed</u> single cells	Spatial Transcriptomics - <u>Visium HD</u>	Spatial Transcriptomics - Visium HD - v2 previous kits: Spatial Transcriptomics - Visium HD	SI-TS-**